

## NEFA Standard Solution

---

### Intended Use

The Wako NEFA Standard Solution is designed to be used to establish the points of reference in the measurement of free fatty acids in serum or plasma, when used with Wako NEFA-HR (2) reagents.

### Reagents and Ingredients

NEFA Standard Solution 4 × 10 mL

Contains : Oleic acid  
Stabilizer

### Determination of assigned values

Assigned values: See bottle label

Determination: Traceable to Wako Internal Standard

### Storage Conditions

Store at 2-10°C

### Calibrator preparation

The standard is supplied in a liquid form and is ready to use.

### Warning and precautions

- \* For In Vitro Diagnostic Use.
- \* Not to be used internally in humans or animals.
- \* This product contains stabilizers. Sodium azide may react with copper or lead plumbing to form explosive compounds. Even though this product contains minute quantities of sodium azide, drains should be flushed well with a large amount of water when discarding the solution.
- \* Do not use this calibrator for any other purposes than described herein.
- \* If the standard comes in contact with the mouth, eye, or skin wash immediately with a large amount of water. Consult a physician if necessary.
- \* When discarding this product, dispose of it according to local or federal regulations.
- \* UDI (Unique Device Identification) is on the box. Keep box until reagent is finished.

### Ordering Information

Code No.	Product	Package
276-76491	NEFA Standard Solution	4 × 10 mL

Revision date ; Apr. 1, 2018

276-76491F

Manufactured by  
**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku Osaka 540-8605, Japan  
Tel : +81-6-6203-3749  
Fax : +81-6-6203-1917  
www.wako-chem.co.jp

**Wako**

Distributed by  
**FUJIFILM Wako Diagnostics U.S.A. Corporation**

1025 Terra Bella Ave., Mountain View, CA 94043 U.S.A.  
Tel : 877-714-1924  
Fax : 804-271-0449  
www.wakodiagnostics.com

18.03.09K02